


Rabbit anti-DiMethyl-Histone H4-K20

Polyclonal Antibody

Catalog Number: R16050P

Immunohistochemistry

	Immunohistochemistry of paraffin-embedded human thyroid cancer tissue using H4K20me2 antibody at dilution of 1:200 (x400 lens)
---	--

Storage

This antibody is shipped at 4 °C. This product is stable for 12 months from date of receipt when stored at -20 °C to -70 °C. Avoid freeze/thaw cycles.

Hazard/Biohazard

This antibody contains 0.02% sodium azide as preservative. Please handle and dispose the product properly. No known biohazard is associated with this product.

For research use only